

# MEDICAL DRESSING SPECIALIST



Healthy Life Aoki Will Be With You Always



Plastic Knitting Industrial Zone, Longgang City, Wenzhou City, Zhejiang, China

+86-0577-64350999/64350988

+86-0577-0577-64350088

www.aoqimed.com

( 🏫

(C,

e





ZHEJIANG AOKI MEDICAL DRESSING CO., LTD





•	Transparent Dressing	P01
٠	Wound Dressing	P03
٠	Hydrocolloid Dressing	P05
٠	Silicone / Alginate And Foam Dressing	P06
٠	Ostomy Bag / Accessories	P07
٠	Hydrogel Dressing	P09
٠	Face Mask	P10
٠	Alcohol Pad	P11
٠	Medical Bandage / IV Catheter Dressing	P12
٠	Medical Tape And Fixation Roll / Wound Skin Closure Strip - · - ·	P13
٠	Surgical Incision Drape / Children Cooling Patch	P14
	Semi Raw Material Roll	P15









Established in 2009 and located in Chinese youngest newborn city Longgang, Zhejiang Aoki Medical Dressing Co.,Ltd is a professional medical consumable supplier which integrate manufacture, sales and research together. Aoki has over 50,000 square meter of ecological garden plant and set up an independently research and development center in Shanghai. With advanced production machine, sterilization equipment and rigorous process technology, Aoki achieve the whole industrial chain cover from production, sterilization and testing.

As a national high technology enterprise, science and technology of Zhejiang Province. Aoki has developed rapidly in the past years by the abundant strength, high quality products and impeccable service system. Involved in three main area, medical, beauty and daily life, Aoki has more than 1000 products category which approved with CE, FDA and ISO13485. By the efforts of the R&D team, Aoki has several patent of invention and utility models.

Aoki will take full advantage of itself to build a new modern plant including well-equipped facilities, reasonable support facilities and thorough function. In order to achieve the goal of world brand based in China in medical dressing industry, Aoki will keep technology, equipment, service and management innovation to satisfy the need in future development of medical dressings.



# COMPANY PROFILE

#### **Transparent Dressing**









B101 4.4cmx4.4cm



A103 5cmx5.7cm







A105 6cmx7cm

B105 6cmx7cm



A107 7cmx8.5cm









A123 10cmx12cm

B123 10cmx12cm









A205 6cmx7cm

B205 6cmx7cm

A405 6cmx7cm

B405 6cmx7cm

#### **Basic Structure**

Type A: PU film, medical adhesive, release paper(film)



#### Features

Transparent dressing adopts a new type of dressing which is ultra-thin and comfortable. It is a special semipermeable membrane with waterproof, breathable, hypoallergic, transparent and elastic. It is widely used in veins and arterial catheters fixation and devices fixation. It can be also used to cover and protect skin damage and provide wound care. It is a necessary product for skin beauty, family standing medicine boxes, military training exercises, wartime first aid, and factory construction.



### Wound Dressing



#### **Basic Structure**

Type B: PU film, medical adhesive, absorbent pad, release paper(film)



#### Features

Transparent dressing adopts a new type of dressing which is ultra-thin and comfortable. It is a special semipermeable membrane with waterproof, breathable, hypoallergic, transparent and elastic. It is widely used in veins and arterial catheters fixation and devices fixation. It can be also used to cover and protect skin damage and provide wound care. It is a necessary product for skin beauty, family standing medicine boxes, military training exercises, wartime first aid, and factory construction.

Model/Specification	Number	Model/Specification	Number
C105 / 6cm x 7cm	50pcs/box 40boxes/carton	C133 / 6cm x 8cm	50pcs/box 40boxes/carton
C120 / 10cm x 10cm	25pcs/box 40boxes/carton	C134 / 6cm x 12cm	25pcs/box 20boxes/carton
C122 / 10cm x 15cm	25pcs/box 20boxes/carton	C136 / 9cm x 20cm	25pcs/box 20boxes/carton
C124 / 10cm x 20cm	25pcs/box 20boxes/carton	C138 / 10cm x 10cm	25pcs/box 40boxes/carton
C125 / 10cm x 25cm	25pcs/box 12boxes/carton	B220 / 10cm x 10cm	25pcs/box 40boxes/carton
C126 / 10cm x 30cm	25pcs/box 12boxes/carton	B222 / 10cm x 15cm	25pcs/box 20boxes/carton
C127 / 10cm x 35cm	25pcs/box 12boxes/carton	B224 / 10cm x 20cm	25pcs/box 20boxes/carton
C132 / 7cm x 9cm	25pcs/box 40boxes/carton	B225 / 10cm x 25cm	25pcs/box 12boxes/carton

\*Can Be Customized According To Customer Requirements













C105 6cmx7cm









C125 10cmx25cm



C132 7cmx9cm











C134 6cmx12cm



C136 9cmx20cm



B220 10cmx10cm



C138 10cmx10cm

B222 10cmx15cm





B124 10cmx20cm

B225 10cmx25cm



B205 6cmx7cm





### Hydrocolloid Dressing

#### Silicone Foam Dressing





#### Characteristics

- Good absorbency for exudation from the wound.
- Keep the wound moist and accelerate the wound healing; reduce the pain and avoid the secondary damage to the wound.
- Waterproof, permeable and prevent the wound from bacteria outside.
- Excellent elasticity and self-adhesive, easy for application.
- Change the dressing when it turns white.
- Easy for peeling off and it will not stick to the wound and keep patients comfortable.



20mmx60mm With border/Extra thin



44mmx69mm With border/Extra thin



82mmx125mm With border



28mmx69mm

With border/Extra thin

50mmx50mm With border/Extra thin



150mmx180mm With border



With border



100mmx100mm

With border/Extra thin

41mmx41mm

With border/Extra thin

150mmx180mm



41mmx62mm

With border/Extra thin

150mmx150mm With border/Extra thin



Acne patch



#### **Characteristics**

Silicone gel application is a soft, elastic application, mainly made of silicon gel, with a certain water vapor permeability and adhesion. Easy to use, and can be repeatedly used. The product can effectively improve and fade the existing scars, prevent the formation of new scars, and help patients regain confidence and beauty.



Shape	Specification	Shape	Specification
	7.5cm x 7.5cm	6	17cm x 18cm
	10cm x 10cm		
	15cm x 15cm		18cm x 18cm

### Alginate Wound Dressing, Foam Dressing

Product	Specification
Alginate Wound	5cm x 5cm
Dressing	10cm x 10cm
Foam Dressing	5cm x 5cm
Foam Dressing	10cm x 10cm

### **Ostomy Bag**










#### **Features**

The hydrocolloid skin barrier has the performance of anti-maceration, low allergy, and rarely warping and leakage. It has excellent exudation absorbent and long sticking ability (up to 7days). It is soft comfortable and convenient to use. The material of stoma bag possesses high obstructing ability, it is softcomfortable and covert and safe.Eliminate odor with active carbon, discharge intestinal gas easily.Humanized design can meet needs of different stoma.



### **Accessories**

Name	Model	Applicable	Name	Model	Applicable
Clamp	81001	Used for sealing the bag	Skin Clean Sheet	81006	Cleaning the stoma wounds arounded
Belt	81002	Strength the bag	Skin Protective Film	81007	Protect the skin from the excrement maceration
Scissors	81003	Cutting the flange in an approriate size	Leakage Proof Strips	81008	Prevent the stoma from leakage
Leakage Proof Cream	81004	Prevent the stoma from leakage	Tweezers	81009	Avoid the touch of the stoma and hand
Ostomy Powder	81005	Absorb the water and keep the skin dried	Stoma Measured Scale	81010	Measure the size in an approriat size



M124





M424



M224

M524

M324

M624

M1



Accessories





### Hydrogel Dressing





Maintain moistness and reduce pain; accelerate wound healing and reduce scars; autolysis and debridement, absorb exudate; Lose stickiness when getting wet, painless when changing; non-toxic, non-irritating, non-allergic; sterile, antibacterial; multi-directional extend, complete and smooth; comfortable and breathable, avoid infiltration; transparent and waterproof, easy to observe.



Shape	Size(mm)	PU Film(mm)	Application
	20 x 20	36 x 36	Laparoscopic surgery, small wounds
	30 x 30	50 x 50	Laparoscopic surgery, small wounds
	50 x 50	70 x 70	Small and medium-sized wounds, surgical wounds, superficial burns and scalds (first-degree and second-degree), chronic wounds
	100 x 100	130 x 130	superficial burns and scalds (first-degree and second-degree), chronic wounds
	25 x 100	50 x 125	Small and medium-sized wounds, surgical wounds, superficial burns and scalds (first-degree and second- degree), chronic wounds
	30 x 140	60 x 170	Caesarean section surgery, surgical wound
	30 x 170	60 x 200	Caesarean section surgery, surgical wound
	10 x 80	25 x 90	Circumcision Surgery
	12.5 x 110	27.5 x 120	Circumcision Surgery

### Face Mask





Non-woven outer layer + Melt-blown fabric filter layer +

Skin-friendly non woven inner layers. Breathable and

Comfortable elastic earloop, extra-soft earloops eliminate pressure to ears. You feel comfortable with

With BFE>98%, perfect when traveling through congested airports, bus terminals, malls, and busy city

- 3-PLY FILTRATION SYSTEM

- COMFORTABLE EARLOOP

periods of time usage.

- STRONG SHIELD

streets.

comfortable, low breathing resistance.



Three fold design The bridge of the nose is bent arbitrarily to support the facial space



Three layers of filtration layer by layer purification



Comfortable and breathable

Product Name	Medical surgical Mask
Model	Non-sterile ear-loop
Specification 17.5cm x 9.5cm	
Number	10pcs/pouch

Product Name	Disposable Face Mask
Model	ear-loop
Specification	17.5cm x 9.5cm
Number	50pcs/box or 20pcs/pouch

### 9



### Alcohol Pad

From Fire



Tableware Disinfection

# Disinfection


### **Features:**

- Suitable for disinfection of the human skin or objects.
- Using Multiple EDI purified water reverse osmosis filtration technology: no impurities.

Disinfection

Sensitive Parts

- Quickly volatilize after wiping, safety in use.
- No artificial flavors, fluorescent agents, preservatives added.
- Spunlace non-woven fabric material: delicate and soft, no remained fibers, not easy to fluff.
- Individual packaging: good leakproofness, easy to carry.



### **Medical Bandage**



Rolls/dozen Dozens/carton

Medical bandage Medical bandage Specification Material Rolls/dozen Dozens/carton Material W

PBT plain cloth

5.0cm x 450cm	12	60
7.5cm x 450cm	12	40
10.0cm x 450cm	12	30
15.0cm x 450cm	12	20

hite wrinkles cloth	, elastic		
	5.0cm x 450cm	12	60
	7.5cm x 450cm	12	40
	10.0cm x 450cm	12	30
	15.0cm x 450cm	12	20

Specification

#### Cotton Crene Bandage

Solion Grepe Bandage					
Council Council	4.8cm x 600cm	12	60		
	6.0cm x 600cm	12	40		
	8.0cm x 600cm	12	30		
$\sim$	10.0cm x 600cm	12	20		

#### Elastic Crepe Bandages

5.0cm x 450cm	12	60
7.5cm x 450cm	12	40
10.0cm x 450cm	12	30
15.0cm x 450cm	12	20


#### Self-adhesive elastic bandage

	Specification	Rolls/box	Boxes/carton
	2.5cm x 450cm	24	24
	5.0cm x 450cm	12	24
	7.5cm x 450cm	8	24
	10.0cm x 450cm	6	24

#### AOKI kinesiology tape

a large app					
	Specification	Rolls/box	Boxes/carton		
⇒	5.0cm x 500cm	12	24		
	7.5cm x 500cm	12	16		

#### IV Catheter Dressing





### Medical Tape And Fixation Roll

Specification

5cm x 10m 10cm x 10m

15cm x 10m

3.5cm x 7cm

3.5cm x 7cm

3.5cm x 7cm

3.5cm x 7cm

Rolls/box

1

1

1

200

40

20

8

Specification Pcs/pouch Pouches/box Boxes/carton

1

5

10

25

Boxes/carton

100

50

50

50

50

50

50

Transparent fixing roll

#### Non-woven fixing roll

Specification	Rolls/box	Boxes/carton
5cm x 10m	1	100
10cm x 10m	1	50
15cm x 10m	1	50

#### Medical tape

Material	Specification	Rolls/box	Boxes/carton
----------	---------------	-----------	--------------

#### Non woven tape

1.25cm x 910cm	24	48
2.5cm x 910cm	12	48
5cm x 910cm	6	48

#### Silk tape

	1.25cm x 500cm	24	30
	2.5cm x 500cm	12	30
	5cm x 500cm	6	30

#### PE Tape

	1.25cm x 500cm	24	48
1000	2.5cm x 500cm	12	48
	5cm x 500cm	6	48

_				
_	PE membrane			
		Acm x Tom	1	

Medical transfusion dressing

Material

Hot rolled non woven

4cm x 7cm	1	100	100
4cm x 7cm	5	20	100
4cm x 8cm	1	100	100
4cm x 8cm	5	20	100

## Surgical Incision Drape



Specification	Pcs/box	Boxes/carton
20cm x 30cm	20	20
30cm x 40cm	20	20
30cm x 45cm	20	20
40cm x 50cm	20	20
45cm x 45cm	20	20



### Children Cooling Patch



Specification	Pcs/pouch	Pouches/box	Bxs/ctn
12cm x 5cm	1	10	40
12cm x 5cm	1	3	400
12cm x 5cm	2	3	400
12cm x 5cm	1	3	400
12cm x 5cm	2	3	400


### Plater patch Base Material

### Wound Skin Closure Strip

Specification	Application	pcs/box	boxes/carton
1.2cm x 10cm	Wound of 5-10cm long	50	40
0.6cm x 10cm	Wound of 5-10cm long	50	40
0.6cm x 3.8cm	Wound of 5cm long	50	40
0.6cm x 7.5cm	Wound of 5cm long	50	40
0.3cm x 7.5cm	Wound of 5cm long	50	40

	Epidermal cells
	Suturing intracutaneous cells
	Suturing subcutaneous cells

subcutaneous cells, make the skin flat with wound skin closures strip.



13



#### Semi Raw Material Roll

Varies materials, for example: Normal type Non-woven Fabric,waterproof non-woven, PU Film and PE Film, PVC Foam ect. with special medical adhesive, low-allergic, strong adhesive and stickness and peel off stength, all of them can be adjusted. Advanced coating technology, with square spot coating, strips coating which can make it more breathable with lower allergic, with leading coating technical in Domestic. All specifications are appliable according to the customer's request.







#### Double CCK paper complex PE adhesive film PU complex PE adhesive film PU adhesive film PE film as backing material. PU film is It weights 140g per square Polyethylene(PE), good 0.025-0.03mm thick, with good ductility ductility, thickness is 0.018-0.02 meter, and peeling strength can and air permeability. Medical pressure be adjusted according to each mm, Medical pressure sensitive sensitive adhesive is 30-32g/m<sup>2</sup>. Colored customer specific request. adhesive is 30-35g/m<sup>2</sup> edge can be sticking automatically while coating. Applied to produce transparent Applied to produce IV catheter dressing, IV catheter fixing fixing dressings, surgical drape, Applied to produce transparent dressing, IV catheter fixing dressing, surgical dressing, surgical drape, and wound plaster base material, drape, wound plaster and band-aids. wound plaster and band-aids. band-aid. PU complex PET adhesive film Hot rolled adhesive non woven Spunlace adhesive non woven (with S cut, easy to be peeled off) spun lance non woven as base Spunlace non woven as base PET as backing material, high material: 45-50g/m<sup>2</sup>, Medical material: 45-50g/m<sup>2</sup>. Medical brightness, good stiffness. Thickness of PU film is 0.025-0.03mm, with good pressure sensitive adhesive: 30pressure sensitive adhesive: 30ductility and permeability, Medical 35g/m<sup>2</sup>, Glassine paper as release 35g/m<sup>2</sup>, Glassine paper as release paper: 100g/m<sup>2</sup>. pressure sensitive adhesive is 30paper: 80-100g/m<sup>2</sup>. 32g/m<sup>2</sup> Applied to produce all kinds of Applied to produce all kinds of Applied to produce transparent dressing, wound dressings, transfusion wound dressings, transfusion IV catheter fixing dressing, surgical drape, wound plaster and band-aids. dressing, wound plaster, heart dressing, wound plaster, heart electrode, etc. electrode, etc. Cuticolor PU acrylic PU film (cartoon-printed) Cartoon PU acrylic complex PET coating film Cuticolor PU film is 0.05mm PET film as backing material, with PU film is 0.05mm thick , with thick, with good ductility and air good stiffness and brightness. PU film good ductility and air is 0.02-0.025mm thick, with good permeability. Printing on PU permeability. Printing on PU ductility and high air permeability. PU membrane. Medical pressure membrane. Medical pressure film is printed with cartoon pictures. sensitive adhesive is 30-35g/m<sup>2</sup>. sensitive adhesive is 30-35g/m<sup>2</sup>. Medical pressure sensitive adhesive is 30-32g/m<sup>2</sup>. Mainly applicable for producing

Mainly applicable for producing new

type band-aids.

new type band-aids.

Mainly applicable for producing new type band-aids.